

# Identification of China's first case of mixed infection of *Chikungunya* virus and *Zika* virus

## MGI Microbiological Detection Total Package

### Highlights

#### Up-to-date Microbial Database

- A comprehensive database of 27, 830 microbial genomics-enabling massive screening at once

#### Proprietary Platforms

- Fully-automated sample preparation system, high-throughput sequencing platform and various supported reagent kits

#### Streamlined Workflow

- Complete the whole process from sample to report with one instrument in 24 hours

#### Flexible Package

- Users can flexibly choose hardware and matching reagent consumables according to their needs

## Introduction

Arboviruses are one of the most common infectious pathogens in the world, among which Dengue virus, *Chikungunya* virus and *Zika* virus can be transmitted through the common host, Aedes. From the analysis of transmission routes, humans are very likely to be infected with multiple arboviruses simultaneously through mosquito bites, which increases the difficulty of clinical diagnosis and treatment.

A passenger departed from Manila, Philippines, to Guangzhou on January 5, 2019, had been having a fever and coughing for 2 days when entering at Guangzhou Baiyun Airport. Epidemiological investigation and blood samples were taken immediately.

The laboratory staff first conducted fluorescence RT-PCR detection of various arboviruses, and the results showed that the nucleic acid amplification of *Chikungunya* virus was positive (Ct value was about 20), while *Zika* virus was positive (Ct value was about 26). It was preliminarily identified that this passenger had mixed infection of *Chikungunya* virus and *Zika* virus. To further confirm the mixed infection results, MGI Microbiological Detection Total Package was used to sequence the virus metagenomic sequence (FIG. 1).




Figure 1. Virus metagenomic sequencing and analysis

## Experimental Method

### Library preparation

The samples were inoculated with Vero cells for virus isolation and culture, and the cells showed obvious lesions after 72 hours. After the nucleic acid was extracted from cell supernatant, the library was prepared using MGIEasy RNA library prep set (16 RXN, Cat. No. 1000006383) and the host rRNA was removed using MGIEasy rRNA depletion kit (32 RXN, Cat. No. 1000005953). The automated sample preparation system MGISP-100 (MGI, Cat. No. 900-000051-00) can also be used to prepare library automatically.

### Sequencing

After each sample library was built, multiple sample libraries with different label sequences were mixed together and paired-end sequencing of 2×50 bp was performed on the DNBSEQ-G50RS\* (Cat. No. 900-000354-00) sequencer. In addition, MGI also provides higher throughput sequencing platforms, DNBSEQ-G400\* (Cat. No. 900-000170-00) and DNBSEQ-T7\* (Cat. No. 900-000128-00), as other options for fast and flexible sequencing packages for pathogen detection.

### Data analysis

MGI has developed a database containing the genetic information of 27, 830 microbes, including bacteria, fungi, viruses, and parasites, and the software named as PFI (Pathogen Fast Identification, Cat. No. 510-000164-00). The integrated system can quickly generate reliable analyses of microbial genome information and the identification report automatically.

## Results

### Verification of mixed infection of *Chikungunya* and *Zika*

The results of viral metagenomic sequencing showed that partial reads were related to *Chikungunya* number of reads was 114,226,369. A small number of reads were related to *Zika*; and the number of reads was 3,320. The results confirmed that the passenger was infected with both *Chikungunya* virus (Table 1) *Zika*

Table 1. PFI software reports the results of virus metagenomic sequencing

NO.	Species	Reads Number	Relative Abundance(%)
1	<i>Chikungunya virus</i>	114,226,369	99.657
2	<i>Barmah Forest virus</i>	282,111	0.246
3	<i>Ross River virus</i>	32,645	0.028
4	<i>Zika virus</i>	3,320	0.003
5	<i>Western quine encephalitis virus</i>	2,052	0.002

## Summary

In this application guide, we introduce the application cases of MGI Microbiological Detection Total Package in the identification of mixed infections. This method combines the library preparation by MGIEasy RNA library prep and the DNBSEQ-G50RS\* sequencer. It also uses the Pathogen Fast Identification software (PFI) for rapid identification of multiple unknown viruses. The whole process (from samples to reports) can be completed in less than 24 hours.

Based on the proprietary high-throughput sequencing platform and rapid identification system of pathogen infection, MGI has developed a microbiological detection total package, which can achieve fast, accurate and comprehensive pathogen screening for clinical diagnoses (FIG. 2). Moreover, MGI provides a variety of hardware and compatible reagent kits for the system to support an extensive range of pathogen testing.




Figure 2. MGI Microbiological Detection Total Package

# MGI

## Microbiological Detection Total Package

### ■ Ordering information

#### Equipment

Product Name	Catalog No.	Configuration
DNA Sequencing Library Preparation System MGISP-100	900-000051-00	RUO
Genetic Sequencer DNBSEQ-G50RS*	900-000354-00	RUO
Genetic Sequencer DNBSEQ-G400RS*	900-000170-00	RUO
Genetic Sequencer DNBSEQ-T7RS*	900-000128-00	RUO

#### Reagents

Product Name	Catalog No.	Configuration
MGIEasy RNA Library Prep Set	1000006383	16 RXN, RUO
MGIEasy RNA Library Prep Set	1000006384	96 RXN, RUO
MGIEasy rRNA Depletion Kit	1000005253	32 RXN, RUO

#### Analysis

Product Name	Catalog No.	Configuration
Fast identification of pathogenic infections (PFI, with server)	510-000164-00	Server + software, RUO
Fast identification of pathogenic infections (PFI)	057-000060-00	software, RUO

MGI Tech Co.,Ltd | Building 11, Beishan Industrial Zone, Yantian District, Shenzhen, CHINA, 518083  
 en.mgi-tech.com | MGI-service@mgi-tech.com

The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

\*Unless otherwise informed, StandardMPS and CoolMPS sequencing reagents, and sequencers for use with such reagents are not available in Germany, USA, Spain, UK, Hong Kong, Sweden, Belgium, Italy, Finland, Czech Republic, Switzerland and Portugal.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

